

# TD-4268

## BLOOD GLUCOSE MONITORING SYSTEM



±15mg/dL & ±15%

### FEATURES

- Strip eject function
- Large LCD with back light
- Display both glucose level and HCT level
- Acoustic mode

5 Sec, 1.0µL

Display Screen



Strip-Ejection



### SPECIFICATIONS

Enzyme Type	GDH
Sample Size	1.0 µL
Reaction Time	5 seconds
Measurement Range	10 ~ 600mg/dL (0.55 ~ 33.3mmol/L)
Hematocrit Range	0% - 70%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100mg/dL (5.55 mmol/L); CV < 5% at 100mg/dL (5.55 mmol/L)
Accuracy	±15mg/dL if < 100mg/dL : ±15% if ≥ 100mg/dL
Ketone Warning	YES
Communication	RS232
Power Source	AAA x 2
Memory Capacity	1000 sets
Day Average	7, 14, 21, 28, 60, 90 days
Daily Alarm	4 daily alarms
Dimension	95 (L) x 55 (W) x 20 (H) mm
Weight	76.2g (Without Battery)
Operating Condition	+10°C(50°F) ~ +40°C(104°F), below 85% R.H.
Storage Condition	-20°C(-4°F) ~ +60°C(140°F) (Meter) ; 2°C(35.6°F)~32°C(89.6°F) (Strip)
Units Per Carton	30
Carton Size	405 (L) x 390 (W) x 365 (H) mm
Cross Weight Per Carton	6.8 kg



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan  
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288