


TD-2555

WEIGHT SCALE

Mb^H
Mobile Health



 **Bluetooth®**



Body mass index function



Age, height, and gender for more precise algorithm



Large LCD with backlight



Automatically power-off function



iOS

Android



Healthy Check



iOS

Android



Healthy Check Pro

SPECIFICATIONS

BMI (Body Mass Index)	Yes
Max. capacity	250Kg
Accuracy	4.0kg ~ 60.0kg ± 0.3kg, 60.1kg ~ 250.0 ± 0.5%
Graduation	0.1Kg
Age Range	10 to 80
Height Range	100 to 220 cm
User	5
Communication	Bluetooth smart
Memory	135 sets
Backlight	LCD Backlight
Power Source	4 x AAA
Dimension	340(L) x 400 (W) x 26 mm (H) mm
Unit per carton	5
Gross Weight per Carton	20kg
Carton Size	48(L) x 29.5 (W) x 43.5 mm (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288