



TD-9030

USB cable

Transfer Measuring Data from Meter to Personal Computer
Support for Window XP or above

Bluetooth Dongle

Transfer Measuring Data from Meter to Smart Phone
Support for iOS and Android platform



SPECIFICATIONS

Bluetooth version	Bluetooth 4.0
Data Interface	Strip Port
Frequency	2.4GHz
Power source	NA (Power by Meter)
Dimension	43.2(L)/15(W)/5(H)mm
Weight	3.56 g
Operating condition	0°C ~ 50 °C (15% ~ 93%RH)
Storage/Transportation condition	-20 °C ~ 60 °C (10% ~ 95%RH)
Units per Box	35
Net Weight per Box	266g
Gross Weight per Box	336g
Box Size	178(L) / 125(W) / 68(H)mm
Applicable models	Blood Glucose / Multi-Functional TD-4141 / TD-4183 / TD-4283 / TD-4206A / TD-4279C

Optionl

Wire transmission	USB
Data Interface	Strip Port



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288